

# Alcohol and Drug Solutions App



Test where and when you need



Easy to input donor information and acknowledgment



Auto Bluetooth Alcolizer/Druglizer results to App



Email results and sync with AlcoCONNECT™



No paper forms, secure electronic platform ensures privacy



Option available for urine test result collection

The screenshot shows a tablet displaying the 'OnSite Testing' app interface. The title is 'OnSite Testing Drug & Alcohol Onsite Screening Form'. There are two main columns of information. The left column is 'Donor Information' and the right is 'Requesting Authority'. Below these are fields for 'Donor Signature Date' and 'Donor's Signature'. At the bottom, there is a line for 'sign above the line'.

Donor Information	Requesting Authority
Name *	Company
Joe Bloggs	Alcolizer HQ
Job *	Collection Site
Leading Hand	Perth
Date Of Birth *	Site Address
01/02/1990	36 Munford Place
Type *	City
Employee	BALCATTA
Reason For Test *	State
Random	WA
Verified By *	Post Code
Two forms of ID	6021
Donor Signature Date	Site Contact
03/03/2022	Roger
Donor's Signature *	
sign above the line	



**OnSite Testing**  
alcohol & drug solutions

Alcolizer Technology  
1300 789 908  
sales@alcolizer.com

The OnSite Testing App is the ideal tool for Safety and Business Managers seeking real-time, mobile enhanced, drug and alcohol test results.

## Why the OnSite Testing App?

Using electronic devices to collect alcohol and drug (AOD) test results via an App means faster testing, secure data and ability to email results and sync with AlcoCONNECT™ Data Management. When combined with AlcoCONNECT™, The Complete Solution, AOD test results can be securely stored in the cloud and at your fingertips when you need.



### RELEVANT, RAPID AND MOBILE RESULTS

Efficiently conduct your AOD testing program where and when you need ensuring your workplace is protected from risks associated with alcohol and drug use. Collect alcohol and drug test results from Druglizer LE5 and Alcolizer LE5 easily using Bluetooth. Minimise risk of data entry error and improve data security. Compliant to AS/NZS 4760:2019 for collection of alcohol and drug test results to improve your testing program.



### EMAIL TEST REPORT OR SYNC WITH ALCOCONNECT™

Results management is as easy as emailing the test report to your nominated address or for The Complete Solution by syncing results to AlcoCONNECT™.



### THE COMPLETE SOLUTION WITH ALCOCONNECT™

When the OnSite Testing App syncs with AlcoCONNECT™ – drug and alcohol test data for all sites is saved in one secure cloud location. AlcoCONNECT™ provides dashboards to display analysed test data and drill down into dashboard data for instant access to results. Dashboards help to identify hot spots in your business and target resources to maximise testing spend.



### SYNC RESULTS AUTOMATICALLY OR ON DEMAND

When internet connectivity is not available, the OnSite Testing App stores results until internet connectivity is restored. Test results will appear in AlcoCONNECT™ when the testing session is completed, signed off and synced manually. Encrypted communication protocol ensures a secure and robust mobile solution with no duplicate records.



### USE APP WITH URINE DEVICE

The OnSite Testing App can be used with a urine drug screening device if required. This provides you with flexibility to use oral fluid or urine to suit your drug testing needs or policy. You can choose from the drop down list the urine device approved for use in the app. Manual entering of the drug test result into the app is required.



### ANDROID AND IOS COMPATIBLE

The OnSite Testing, Alcohol and Drug Solution App is available for use with Android and iOS operating systems. The download links and instructions to download are available directly from Alcolizer Technology. Contact your local sales representative for more information.

## Overview of the App

The App receives screening drug and alcohol test records from a paired Bluetooth or Bluetooth Low Energy (BLE) enabled Alcolizer and/or Druglizer handheld device.

Non-Druglizer, urine and oral fluid drug screening devices can be chosen from the dropdown list of approved devices and test results can be manually entered into the App.

Donor details are entered manually into the app in line with Australian Standard AS/NZS 4760:2019 requirements and declarations are acknowledged via electronic signing within the app.

These test records can be synced with AlcoCONNECT™ or emailed as PDF documents directly from the app.




Image above: Donor information is entered into the app




Image above: The OnSite Testing App and AlcoCONNECT™. Results are received via Bluetooth from Alcolizer or Druglizer devices and can be synced with AlcoCONNECT™.





Image above: Test results from Alcolizer or Druglizer devices are received automatically when devices are paired with the app. Image to the right depicts a DOB selection.



# AlcoCONNECT

data management

Live Drug & Alcohol Test Data. AlcoCONNECT™ combines Alcolizer’s innovative testing technology with state-of-the-art business solutions.

- ✓ Simple dashboards with real-time data analysis
- ✓ Drill-down to a single test result with photo evidence\*
- ✓ Monitor compliance in daily testing programs
- ✓ Secure cloud infrastructure
- ✓ Target your activity to maximise return on investment
- ✓ Real-time analysed alcohol and drug test data from across your business.

\* When camera option is available and activated. Camera is not applicable to handheld devices/OnSite Testing app.

GPS map location is not available for the OnSite Testing app.



Contact us for all your training and education requirements

1300 789 908  
 sales@alcolizer.com  
 www.alcolizer.com

V 700.01.250322

Alcolizer  
 technology

TRUSTED ALWAYS